



NURIN^S

A Nurin user guide ————— *version 21*
With introductions to Nurin Nurin/Nurin Yusof/Tengku Nurin/Nurin Girl Child/Nurin Exotic Asian Girl/Nurin Imaginary Friend

CONTENTS

1	Introduction	3
	1.1 Nurin Nurin	3
	1.2 Nurin Yusof	3
	1.3 Tengku Nurin	3
	1.4 Nurin Girl Child	3
	1.5 Nurin Exotic Asian Girl	3
	1.6 Nurin Imaginary Friend	3
2	Product Diagram	4
3	How to Use	5
	3.1 General Functions	6
	3.2 Assembly	7
	3.3 Maintenance	8
	3.4 Changing Modes	9
	3.4 Additional Accessories	10
4	Personalisation	11
5	FAQ	12
	5.1 General FAQ	13
	5.2 Troubleshooting Nurin	14
6	Contact	16

(5.1) Troubleshooting Nurin

(1) Broken Tear Gland(s) (levels):

Level 1 - Teary Eyes; water accumulating in the corner of eye(s)
Level 2 - Single Tear; single water droplet falling from single tear gland
Level 3 - Crying; repetitive spillage of water droplets from tear glands
Level 4 - Sobbing; crying with choking sounds due to lack of breathing
Level 5 - Heaving; sobbing with wailing sounds and physical movement
Level 6 - Silent Crying/Sobbing/Heaving; sound suppressed due to significant physical and/or emotional pain

Level 7 - Silent Crying/Sobbing/Heaving in bath or shower; naked
Level 8 - Silent Crying/Sobbing/Heaving in room corner or cold toilet; comfort in location
Level 10 - Crying/Sobbing/Heaving in bed, before sleep or to fall asleep
Level 11 - Silent Crying/Sobbing/Heaving to Nurin's reflection; red eyes reflecting desperation and sincerity of Nurin as the self
Level 12 - Crying/Sobbing/Heaving whilst consciously talking to Nurin's reflection

Level 13 - (Silent) Crying/Sobbing/Heaving in public; self control displaced
Level 14 - Crying to Aminah at The Coffee Bean & Tea Leaf, Bangsar Village II, outside seat by the window

To Fix:

Step 1 - Remove Nurin from situation (location, family, friend(s), school, etc.)
Step 2 - Give Nurin Topsy
Step 3 - Leave Nurin
Step 4 - Give Nurin phone to contact Aminah or for Nurin to contact Nurin
Step 5 - Place Nurin in front of mirror for conscious rebooting
Step 6 - Undress Nurin, then place Nurin in hot bath for 2 hours
Step 7 - Remove Nurin from bath, pat dry
Step 8 - Place Nurin in bed and leave to sleep

Rebooting will commence once awake and will take a minimum of 24 hours to complete.
With each increasing level, more time is needed for rebooting and may require the steps to be repeated and altered accordingly.

(2) Overheating of Adrenal Glands (Level):

Level 1 - Silence with neutral facial expression

Level 2 - Silence with smiling facial expression

Level 3 - Silence with heavy breathing; increase in heart rate, arterial tension and respiration

Level 4 - Conscious Violent Acts; throwing of personal objects, pinching and scratching of self and pulling of hair from roots

Level 5 - Unconscious Angry-Helpless Tears

To Fix Level 1 — Level 3:

Step 1 - Remove Nurin from situation (location, family, fiend(s), school, etc.)

Step 2 - Place Nurin in a room alone and allow Nurin to walk in circles whilst talking to herself

Step 3 - Give Nurin phone to contact Aminah or for Nurin to contact Nurin

Step 4 - Place Nurin in front of mirror for conscious rebooting

Step 5 - Sit Nurin down on the floor in the corner of the room to steady breathing, then leave to cool

To Fix Level 4:

Step 6 - Complete *step 1 to step 5*, repeat if necessary and if problem still arises, proceed to *Step 7*

Step 7 - To halt violent acts, remove personal belongings from Nurin's hands

Step 8 - If Nurin continues to act violently, restrain Nurin and replace objects with soft objects

Step 9 - Allow Nurin to complete violent acts once the objects are replaced

Step 10 - If Nurin begins to pinch and scratch Nurin, restrain Nurin and cut Nurin's nails

Step 11 - If Nurin begins to pull her hair, cover hair temporarily and give Nurin soft objects to hold

Step 12 - Once soft object is Nurin's possession, remove hair cover and repeat *Step 5*

To Fix Level 5:

Step 13 - Complete *step 1 to step 5*, repeat if necessary and if problem still arises, proceed to *Step 14*

Step 14 - Identify Tear Level by referring to *Broken Tear Gland(s)* troubleshooting guide and proceed to take action accordingly

Cool down will commence once breathing is steadied. This could take a minimum of 15 minutes, but could exceed over 24 hours. What should be noted is that Nurin should be left in silence throughout the cool down process.